


# TRANSFERRING PATIENT TRENDELENBURG

# 31209T



## CRANK HANDLE

Equipped with excellent crank driven enable full adjusted product solution



## UPHOLSTERY

Provided upholstery with positional wide surface excellent support and comfort for patients treatments.



## BASKET

Equipped excellent stainless steel holder basket.




## CASTORS

Giving smooth product moving fitted 6 inch double wheels castor dan 5 inch steering wheel.

SPECIFICATION	
Construction	High quality steel sections construction utilizing the latest manufacturing techniques throughout.
Mattress Platform	Manufactured from perforated steel sheet.
Upholstery	Manufactured in line using high quality foam with vinyl.
Hi - Lo	Adjustable manually by crank, elevation 700 - 930mm.
Back-raise	Inclination 0 - 40° - gas spring mechanism.
Side guard	Drop down ABS Side guard.
Trendelenburg	Inclination 0 - 12° - Gas spring mechanism.
Accessories	- O2 Bracket - Stainless steel basket - Stainless steel section infuse pole.
Castors	6" Central lock double wheel castors, and 5" wheel.
Coating	Epoxy powder coating finish.
Max. Load	170 kg.
Net Weight	± 108 kg.
Dimension	(L)2050 x (W)680 x (H)890 - 1120 mm (overall). (L)1915 x (w)600 x (H)700 - 930 mm (on mattress platform)

Specifications subject to change without prior notice.

## DIMENSIONS (mm)



Eff. Date 02-09-2024



Marketing Office :  
Gran Rubina Business Park Lt. 22C. Komp. Rasuna Epicentrum, Jl. HR Rasuna Said - Jakarta 12940  
Ph. +62-21-837 00 555. Fax. +62-21-837 00 355, Email: mak@cbn.net.id/marketing@mak-techno.com  
www.mak-techno.com